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- Do not use sampling device if broken or damaged.
- Do not use test cartridge or sampling device after the expiration date.
- Read results after exactly 20 minutes. Do not interpret results after 25 minutes. Depending on the concentration of the biological threat agent, an earlier, qualitative, visual read-out is possible. Positive results might be already observed after 2 minutes.
- All materials used have to be correctly disposed of.

#### Storage and stability

Test cartridge should be stored cooled or at room temperature (+4°C to +30°C). During storage, frequent changes in temperature should be avoided. Test cartridge must not be frozen. It must remain in foil pouch until use. Sampling device as well as test cartridge must not be used after the expiration date printed on the foil pouch (test cartridge).

#### Special note and limitations



1. A negative test result does not exclude a possible biological threat agent but reduces the possibility of a threat.
2. In case of a positive test result local regulations and laws must be observed.
3. miPROTECT® rapid test provides a preliminary analytical result. Confirmation should be made by an additional method.
4. There is always the possibility that factors such as technical or procedural errors as well as interfering substances in the sample may cause erroneous results.

#### Sample collection and test procedure

Remove the test cartridge from the foil pouch and place it on a dry, clean, flat surface. If required, write the sample identification on the intended field (ID) of the test cartridge. After opening the foil pouch and removing the test cartridge, the test procedure has to be started without further delay.

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#### Sample collection from surfaces or solids (powder)

- ▶ Snap the buffer reservoir until the coloured little pin inside is broken. **(A)**
- ▶ Squeeze the buffer reservoir several times with two fingers to get the buffer to flow into the tube to moisten the swab. **(B)**
- ▶ Remove the swab from the tube and skim the swab against the inside of the tube to remove excess liquid from the swab. Do not touch the swab.
- ▶ Rotate and rub the moistened swab in the sample to collect sample material. **(C)**
- ▶ Place the swab with the collected sample into the tube containing the buffer and reassemble both parts.
- ▶ Rotate the swab several times (10 times) inside the tube to remove the sample material from the swab. **(D)**
- ▶ Test procedure see back

#### Sample collection from a liquid

- ▶ Remove the swab from the tube. Do not touch the swab.
- ▶ Rotate and rub the swab in the liquid sample to collect sample material. **(E)**
- ▶ Place the swab with the collected sample into the tube and reassemble both parts.
- ▶ Snap the buffer reservoir until the coloured little pin inside is broken. **(A)**
- ▶ Squeeze the buffer reservoir several times with two fingers to get the buffer to flow into the tube to moisten the swab. **(B)**
- ▶ Rotate the swab several times (10 times) inside the tube to remove the sample from the swab. **(D)**
- ▶ Test procedure see back

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## Test procedure

Remove the white cap from the bottom of the sampling device.

Add 4 drops of the diluted sample onto the oval sample port (S) of the test cartridge by slightly squeezing the tube of the sampling device. The tube has to be held at a slight angle.

Do not add more than 4 drops onto the test cartridge.

Re-close the sampling device with the white cap.

Interpret test results after exactly 20 minutes. For qualitative visual read-out, tests can be read out earlier depending on the concentration of the biological threat agent. Positive results might be already observed after 2 minutes. Do not interpret test results after 25 minutes.

## Interpretation of results





Read results in the result window after exactly 20 minutes.

**Negative:** The appearance of one color line (control line) next to the letter C in the result window indicates a negative result.

**Positive:** The test is positive if two colored lines appear. One colored line will appear next to the letter C (control line) and another colored line (test line) will appear next to the letter (T).

**Invalid:** The test is invalid if no colored line appears.

**Invalid:** The test is invalid if a test line appears next to the letter (T), only.

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Rapid test for the detection of biological threat agents

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*Please read the following guidance note carefully.*

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## Intended use

miPROTECT® is an easy-to-use rapid test for testing of environmental samples. The test is suitable for a visual read-out. In addition, tests can be measured with a portable reader provided by mipro lab. This reader allows a semi-quantification of the test signals. The test is intended for screening of environmental samples and NOT for any medical use.

## Test principles

miPROTECT® is a rapid test for the detection of biological threat agents in environmental samples. After applying the diluted sample onto the sample port (S) of the test cartridge, the solution migrates via capillary forces through the test cassette. After less than a few minutes, a red line appears in the control zone (C) in the result window. This line confirms that the test is working correctly. Another red line appearing in the test zone (T) shows a positive signal and the presence of the biological threat agent in the sample. A negative test consists of a red line appearing solely in the control zone (C).

## Materials provided

5 test cartridges (separately sealed in a foil pouch)  
5 sampling devices (ready-to-use)  
1 package insert

## Materials required but not provided

Timer  
Disposable gloves  
Appropriate personal protective equipment

## Precautions

- For use by qualified persons only.
- Handle all samples as extremely toxic and as potentially infectious.
- Wear disposable gloves and appropriate personal protective equipment.
- Do not eat, drink or smoke while handling the specimens.
- Test cartridge should remain in the sealed foil pouch until use. Foil pouch should be opened only just before test procedure.
- Do not use test cartridge if foil pouch is damaged or open.
- Test cartridge is for single-use, only. Do not re-use test cartridge.